

TRENDGOLD Dentallote  CE0123	Farbe	Arbeits- temperatur in °C	Schmelz- intervall in °C	Massenanteile in % (x < 1%)										Gewicht in g	Indikationen
				Au	Pd	Pt	Ag	Cu	Zn	In	Ni	Ir	Sonstige		
Lot 0 850	gelb	850	810-860	72,0	x	2,0	10,5	9,2	5,4	-	-	-	-	5	5
Lot 1 800	gelb	800	760-800	72,0	1,0	1,9	8,6	7,1	9,4	-	-	-	-	5	5
Lot 2 750	gelb	750	710-760	73,2	-	1,8	10,6	3,0	11,4	-	-	-	-	5	3,4,5
Lot 1 PF 800	gelb	780	690-800	70,5	-	1,4	6,5	13,0	3,0	5,5	-	x	-	5	3,4,5
Lot N1 PdF 880	gelb	880	820-880	76,0	-	2,9	10,0	6,0	5,0	-	-	x	-	5	2
Lot N2 PdF 700	gelb	700	650-710	72,5	-	x	10,0	3,0	12,0	2,0	-	x	-	5	3,4,5
Lot Uni 1	gelb	870	760-870	56,0	-	-	36,0	-	8,0	-	-	-	-	5	5
Lot Uni 2	gelb	800	690-800	53,0	-	-	35,5	-	11,5	-	-	-	-	5	5
Lot 710	hell- gelb	710	670-720	50,0	1,0	-	27,5	5,1	14,0	2,3	-	-	Fe x	5	3,4,5
Lot Alpha 1 950	gelb	960	930-970	60,5	-	x	35,1	-	1,9	-	-	x	Sn 1,9	5	2
Lot Alpha 2 880	gelb	880	780-890	80,0	-	2,0	11,4	-	6,5	-	-	x	-	5	2
Lot Alpha 3 680	gelb	690	630-710	72,4	-	x	10,2	2,9	12,0	2,0	-	x	-	5	3,4,5
Lot U1 1120	weiß	1140	1030-1140	70,4	9,7	-	17,5	x	1,0	1,0	-	-	-	5	1
Lot N1 1060	gelb	1070	980-1070	80,0	6,0	x	8,5	3,2	1,9	-	-	-	-	5	1
Lot UW 1060	weiß	1060	950-1060	71,9	12,0	-	5,0	8,0	2,0	1,0	-	x	-	5	1
Lot G1 1030	gelb	1060	990-1070	72,0	7,6	-	16,0	1,6	2,8	-	-	-	-	5	1
Lot T1-PdF 1040	gelb	1050	980-1050	79,0	-	3,0	16,6	-	1,3	-	-	x	-	5	1
Lot Bio Plus 1010	gelb	1020	1020-1050	64,0	-	x	34,9	-	-	-	-	x	Rh x	5	1
Lot StG 750	blass- gelb	770	670-770	42,0	3,0	-	24,4	16,5	4,0	7,0	-	x	Mn 3,0	5	6
Lot StG 930	weiß	930	880-930	79,8	-	-	-	-	4,4	-	15,8	-	-	5	6

1=Vor-Brand-Lot

2=Vor-Brand-Lot LFC- Legierungen

3=Nach-Brand-Lot

4=Nach-Brand-Lot LFC-Legierungen

5=Verbindungslot Gusslegierungen

6=Verbindungslot NEM/EM-Legierungen